



FEMTOCHROME[®]
RESEARCH, Inc.

2485 4th Street., Berkeley, CA 94710
e-mail: sales@femtochrome.com;

Ph#: 510-644-1869
http: www.femtochrome.com

FR-103TPM AUTOCORRELATOR



Specifications

- Sensitivity $(P_{av}P_{pk})_{min} : 10^{-4}W^2$
 - Resolution: $\sim 1fs$
 - Scan Range: $> 50 ps$
 - Wavelength Range: 700–2200nm
 - Interferometric
 - Polarization Insensitive (TPC)*
 - Computer Interface (/CDA)
 - Any pulse rep rate $> 4Hz$ (w/CDA)
- * Vertical polarization provides higher R/T ratio for pellicle beamsplitters

The **FR-103TPM** is a compact autocorrelator for Two Photon Microscopy applications. It can be introduced into the beampath of the laser without affecting optical alignment. The NL photosensor of the FR-103TPM can be placed in any location of the optical set-up, in particular at the position of the sample. Its Computer Data Acquisition (/CDA) option provides an interface (USB), for the display and analysis of autocorrelation traces on a Windows PC.

The standard **FR-103TPM** works with any input pulse rep rate $> 500Hz$. [This can be extended to any rep rate $> 4Hz$, by the addition of /SSO and/CDA options.]


DISPERSION-FREE, HIGH RESOLUTION

Dispersion is negligible in the **FR-103TPM**. Using high reflective metallic-coated optics [the only transmissive element is an ultrathin ($<1\mu m$) pellicle beamsplitter], an unprecedented resolution of $\sim 1fs$ is attained.

With its high resolution, the **FR-103TPM** is capable of measuring pulsewidths as low as $\sim 5fs$.

ROTATING PARALLEL (//) MIRROR ASSEMBLY

Rapid scan, periodic optical delay is introduced by means of a parallel (//) mirror assembly.* This unique mechanism results in uniform and error-free delay generation with interferometric resolution.



* Z.A.Yasa and N.M.Amer, Optics Commun., V36, 406 (1981).



NONLINEAR PHOTSENSOR MODULES

The NL detector module selections for the **FR-103TPM** are:

- /700 → 700-1200nm
- /1200 → 1200-2200nm

SLOW SCAN OPERATION

There are two modes of operation of the **FR-103TPM**:

1. Uniform rotation of // mirrors(with a refresh rate typically ~2Hz)
2. Controlled movement such that the // mirrors slow down greatly (4 selectable speeds) over a period when the pulses on the two arms of the Michelson Interferometer set up are overlapping. The // mirrors speed up beyond this range, to return quickly, for a repetition of the cycle.


Whereas the uniformly rotating // mirrors is suitable for typically > 100kHz rep rate lasers, the latter mode renders the unit suitable for 'real-time' autocorrelation for any rep rate > 500Hz. This is particularly useful for kHz amplified lasers.

COMPUTER DATA ACQUISITION (/CDA)

A data acquisition board is installed in the **FR-103TPM/CDA**, provides a USB interface with any PC w/ Windows OS. Its associated software allows traces to be displayed, analyzed [averaged and/or fit with typical pulseshapes (Gaussian and Sech^2)] or saved.

SPECIFICATIONS:

- * Resolution: ~ 1fs
- * Scan Range: > 50ps
- * Sensitivity: $[P_{av}P_{pk}]_{min}=10^{-4}W^2$
- * Wavelength Range: 700-2200nm
- * Interferometric
- * Polarization insensitive*
- * Any rep rate > 4Hz (w/CDA)
- * Computer Interface (/CDA)
- * Vertical polarization generally yields higher signal



2485 4th Street., Berkeley, CA 94710
e-mail: sales@femtochrome.com;

Ph#: 510-644-1869
http: www.femtochrome.com