



Datasheet PRA-PAP

Post Adapter Plate for Photoreceivers



Features	 Compatible with all FEMTO photoreceiver series FWPR, OE, LCA-S and HCA-S M4 and 8-32 UNC threads suitable for standard optical mounting posts High-tensile material Mounting screws included
Applications	 Mounting FEMTO photoreceivers to optical tables/breadboards Expanding the mounting options for all FEMTO photoreceiver series FWPR, OE, LCA-S and HCA-S
Scope of Delivery	PRA-PAP, 4 x M2.5 x 12 mounting screws, datasheet, transport package
Specifications	
General Data	Material AIMg3; glass bead blasted, nickel-plated Weight 38 g (0.084 lb.)
Dimensions	A

FEMTO Messtechnik GmbH Klosterstr. 64 10179 Berlin · Germany Phone: +49 30 280 4711-0 Fax: +49 30 280 4711-11 Email: info@femto.de www.femto.de Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

© by FEMTO Messtechnik GmbH · Printed in Germany


SOPHISTICATED TOOLS FOR SIGNAL RECOVERY





DE-PRA-PAP_R1/MG,JM/15APR2015 Page 1 of 2

Mounting Instructions

PRA-PAP







Choose position.

You can mount the adapter at position A, B or C at the bottom of the photoreceivers of the OE series.


At the photoreceivers of the series FWPR, LCA-S and HCA-S the adapter plate must be mounted to the center position.

Take off screws.

Four screws have to be taken off at the chosen position. Do not take off more than four screws! Do not take off screws at the corners!

Choose thread.

Put the plate in the choosen postion. Turn the plate to position the relevant thread (M4 or 8-32 UNC) in line with the optical axis.



Fasten adapter plate.

Take the four M2.5 x 12 screws, which were provided with this mounting set, and fasten the post adapter plate.

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY