



Instrument Expert Original factory packaging www.dorgean.com

Uniblitz[®] XRS14

14mm Uni-Stable X-ray Shutter

Overview

The Uniblitz XRS14 has been designed specifically for x-ray switching applications. The innovative platinum-iridium (PtIr) blade design allows beam extinction of 90% up to 30 keV x-ray energy (based upon the tenth value extinction of PtIr). The XRS14 is well suited for applications such as x-ray crystallography, medical x-ray imaging, etc.

Uni-stable shutter devices, like the XRS14, require power to hold the blades in the open state.

Need Support? Please visit our website or email us at info@uniblitz.com.

Tel: <u>585-385-5930</u> | Toll-Free: <u>800-828-6972</u> | Fax: 585-385-6004 | 803 Linden Ave. Rochester, NY 14625 Updated 11/18 | Datasheet Version 5.2 | ©2018 Vincent Associates




Key Features

- 14mm aperture
- Uni-stable operation
- Capable of blocking x-ray energy (30 KeV)
- RoHS Compliant
- Transfer time on opening:
 - 20.0 milliseconds
- Total opening time:
 - 25.0 milliseconds
- Configured for the <u>VCM-D1</u>
 <u>Shutter Driver</u>

Product Options

XRS14 2 3 4 5 6	- 7 - 8 Ex: XRS14S2P0-EC-21
Shutter Series:XRS14	 2 Driver Compatibility: S: Use with VCM-D1 (Std.) E: Use with D880C or VED24
 3 Housing: 1: Un-Housed 2: #2 Housing 	 Blade Coating: P: 0.01" thick Ptlr (10% iridium, 90% platinum)
 Electronic Sync: 0: Omitted 1: Included 	 Connector: L: 18" flying leads (Un-housed only) Leave blank for 7-pin Wire Pro connector
 Encapsulated Coil: EC: Included Leave blank if not required 	 8 Mount: (#2 housing only) 21: Zeiss Axiovert 24: Olympus BX/IX 27: Nikon 105: C-Mount (Male) 106: C-Mount (Female) Leave blank if not required

Shutter Timing



XRS14 (w/ VCM-D1 driver and PtIr blades)Time (msec.)O - ADelay time on opening after current applied5.0A - CTransfer time on opening20.0O - CTotal opening time25.0C - EMin. dwell time with min. input pulse20.0

25.0	Total opening time	0 - C
20.0	Min. dwell time with min. input pulse	C - E
5.0	Min. equivalent exp. time	B - F
10.0	Transfer time on closing	E - G
35.0	Total window time	A - G
25.0	Min. exposure time	MET
>25.0	Typical exposure pulse	TEP


Technical Specifications

Coil Resistance	Voltage to Open	Hold Voltage (Nominal) ¹	
12 Ω	+65 VDC	+5 VDC	

¹ Voltage level required across actuator coil when being held in the open position.


² (Continuous/Burst) Continuous frequency rating specified at shutter's minimum exposure pulse. Burst frequency rating specified for four (4) seconds maximum with one (1) minute minimum between bursts.

Series	Weight	Weight	Operating	Max. Opening	Max. Closing	Max. Freq. of	Number of
	(Unhoused)	(Housed)	Temp.	Bounce	Bounce	Operation ²	Shutter Blades
XRS14	2.18 oz (0.06 kg)	10.35 oz (0.29 kg)	0 - 80 °C	15%	5%	2 Hz / 10 Hz	1




The electronic synchronization system provides a feedback signal (through the driver utilized) after the shutter transfers to the open state. The system incorporates an infrared emitting diode, an infrared sensitive detecting transistor, and an interrupting vane. The vane is attached to the shutter so as to block the light path between the emitter and detector in the closed position. When the shutter transfers to the 80% open position, the vane is removed from the infrared light path, allowing the emitter to switch the detector to the active state. For the XRS14, this system uses a similar activation flag attached to the mechanism, which triggers a reflective emitter/detector device. **No connection to the designated synchronization pins when no electronic sync. is selected.**

Technical Drawings - Un-housed XRS14



Technical Drawings - Housed XRS14



NOTES:

- (3) #2-56 MOUNTING HOLES ON A Ø 1.204 BOLT CIRCLE, 120° APART
 (3) #2-56 MOUNTING HOLES ON A Ø 1.578 BOLT CIRCLE, 120° APART
- 2. (3) #4-40 MOUNTING HOLES ON A Ø 2.425 BOLT CIRCLE, 120° APART (UNDER NAME PLATE)
 - ALL DIMENSIONS [MM] INCH



Technical Drawings - XRS14 Connector Layout

