



# gentec-()

## UP52M-300W-QED-D0

Thermal detector for laser power measurement up to 300 W.



## PRODUCT FAMILY KEY FEATURES

## MODULAR CONCEPT

Increase the power capability of your detector: 4 different cooling modules

HIGH PEAK POWER DIFFUSING ABSORBER

Perfect for pulsed beams with high energy density

40 mm thick

HIGH AVERAGE POWER Measure up to 300 W of continuous power

SMART INTERFACE Containing all the calibration data

#### AWARD-WINNING TECHNOLOGY

The UP-QED laser power detectors for extremely high density lasers were recognized among the most innovative photonics technologies for the 2021 Laser Focus World Innovators Awards, as a Gold honoree.



COMPATIBLE STAND

STAND-S-443



### MEASUREMENT CAPABILITIES

Maximum average power (continuous) <sup>1</sup>	300 W
Maximum average power (1 minute) <sup>2</sup>	300 W
Noise equivalent power <sup>3</sup>	15 mW
Spectral range <sup>4</sup>	0.266 - 2.5 μm
Typical rise time <sup>5</sup>	4 s
Power calibration uncertainty <sup>6</sup>	±2.5 %
Repeatability	±0.5 %

1. Minimum cooling flow 1 liters/min, water temperature < 22°C, 1/8 NPT compression fittings for 1/4 inch semi-rigid tube. Contact Gentec-EO for clean deionized water cooling module option.

2. Minimum cooling flow 1 liters/min, water temperature < 22°C, 1/8 NPT compression fittings for 1/4 inch semi-rigid tube. Contact Gentec-EO for clean deionized water cooling module option.

3. Nominal value, actual value depends on electrical noise in the measurement system.

1. For 360 µs pulses. Higher pulse energy possible for long pulses (ms), less for short pulses (ns).

For the calibrated spectral range, see the user manual.
With anticipation.

6. Including linearity with power.

## MEASUREMENT CAPABILITIES (ENERGY MODE)

Maximum measurable energy <sup>1</sup>	1000 J
Noise equivalent energy <sup>2</sup>	0.25 J
Minimum repetition period	9 s
Maximum pulse width	371 ms
Energy calibration uncertainty <sup>3</sup>	±5 %





Nominal value, actual value depends on electrical noise in the measurement system.
When single-shot energy calibration is purchased

DAMAGE THRESHOLDS	
Maximum average power density <sup>1</sup>	100 kW/cm <sup>2</sup>
Maximum energy density <sup>2</sup>	8 J/cm²
1. At 1064 nm, 10 W CW. May vary with wavelength and average power. 2. At 1064 nm, 7 ns, 10 Hz. May vary with wavelength and pulse width.	
PHYSICAL CHARACTERISTICS	
Cooling	Water
Aperture diameter	52 mm
Absorber	QED
Dimensions	89H x 89W x 43D mm
Weight	0.84 kg
ORDERING INFORMATION	
UP52M-300W-QED-D0	203878
UP52M-300W-QED-IDR-D0	205205
UP52M-300W-QED-INT-D0	205198
UP52M-300W-QED-BLU-D0	TBD

Specifications are subject to change without notice. Refer to the user manual for complete specifications.

## **INTERESTED IN THIS PRODUCT?**

GET A QUOTE

Find your local sales representative at gentec-eo.com/contact-us