

Datasheet HSA-X-S-1G4-SI

Ultra High Speed Photoreceiver with Si-PIN Photodiode



The picture shows model HSA-X-S-1G4-SI-FS. The photoreceiver will be delivered without post holder and post.


Si-PIN photodiode
 Bandwidth 10 kHz – 1.4 GHz
 Amplifier transimpedance gain 5.0 × 10³ V/A
 Max. conversion gain 2.55 × 10³ V/W @ 760 nm
 Spectral range 320 – 1000 nm

- Free-space input 1.035"-40 threaded, alternatively 25 mm diameter unthreaded
- Easily convertible to fiber optic input (FC and FSMA) with optionally available screw-on adapters
- UNC 8-32 and M4 tapped holes for mounting on standard posts with metric and imperial thread

Applications

- Spectroscopy
- Ultra-fast pulse and transient measurements
- Optical triggering
- Optical front-end for oscilloscopes and ultra-fast A/D converters

Block Diagram



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E 🏲 T O

HSA-X-S-1G4-SI_R11/TH,JMa/31MAR2022 Page 1 of 7

Ultra High Speed Photoreceiver with Si-PIN Photodiode

Available Versions

HSA-X-S-1G4-SI-FST



Picture shows 1.035"-40 threaded flange with internally threaded coupler ring (outer diameter 30 mm)

1.035"-40 threaded flange for free space applications. Compatible with many optical standard accessories and for use with various types of fiber connector adapters.

Optionally available:

Fiber adapters PRA-FC, PRA-FCA and PRA-FSMA. The coupling efficiency will depend on fiber type. With the relative large 0.4 mm dia. photodiode installed in the HSA-X-S-1G4-SI input coupling is not critical. However, standard SM 9/125 fibers (PC or APC) with low numerical aperture (NA) are recommended for ensuring near 100% coupling efficiency.

HSA-X-S-1G4-SI-FS



Picture shows unthreaded flange with 25 mm diameter

25 mm dia. unthreaded flange for free space applications. Compatible with many optical standard accessories.

HSA-X-S-1G4-SI-FC



Fix/permanent FC fiber connector for high coupling efficiency and excellent conversion gain accuracy.

Related I	Models
-----------	--------

HSPR-X-I-1G4-SI-FST

Si-PIN, \varnothing 0.4 mm, 320 – 1000 nm, inverting output free space input, 1.035"-40 threaded flange

HSPR-X-I-1G4-SI-FS

Si-PIN, \varnothing 0.4 mm, 320 - 1000 nm, inverting output free space input, 25 mm dia. unthreaded flange

HSPR-X-I-1G4-SI-FC

Si-PIN, integrated ball lens, 320 – 1000 nm, inverting output, FC fiber connector (fix/permanent)

HSA-X-2G-IN-FST

InGaAs-PIN, \varnothing 0.1 mm, 900 – 1700 nm free space input, 1.035"-40 threaded flange

HSA-X-2G-IN-FS

InGaAs-PIN, Ø 0.1 mm, 900 − 1700 nm

HSA-X-2G-IN-FC

free space input, 25 mm dia. unthreaded flange InGaAs-PIN, integrated ball lens, 900 – 1700 nm

FC fiber connector (fix/permanent)

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

Е 🍑 Т О

HSA-X-S-1G4-SI_R11/TH,JMa/31MAR2022 Page 2 of 7

Datasheet

HSA-X-S-1G4-SI

Ultra High Speed Photoreceiver with Si-PIN Photodiode

Related Models (continued) HSPR-X-I-2G-IN-FST

InGaAs-PIN,Ø 0.1 mm, 900 − 1700 nm, inverting output

free space input, 1.035"-40 threaded flange

HSPR-X-I-2G-IN-FS InGaAs-PIN, Ø 0.1 mm, 900 − 1700 nm, inverting output

free space input, 25 mm dia. unthreaded flange

HSPR-X-I-2G-IN-FC InGaAs-PIN, integrated ball lens, 900 – 1700 nm, inverting

output, FC fiber connector (fix/permanent)

Available Accessories

PRA-FC PRA-FCA PRA-FSMA







Fiber-adapter with external 1.035"-40 thread (suitable for FST models only).

PS-15-25-L



Power supply Input: 100 - 240 VAC Output: ±15 VDC

Specifications

Test conditions

 $V_S = +15 \text{ V}$, $T_A = 25 \,^{\circ}\text{C}$, output load impedance 50 Ω ,

warm-up 20 minutes (min. 10 minutes recommended)

Gain

Transimpedance gain Conversion gain

 5.0×10^3 V/A (@ output load 50 Ω)

 2.55×10^{3} V/W typ. (@ 760 nm, output load 50 Ω)

Frequency Response

Lower cut-off frequency (-3 dB)

10 kHz

Upper cut-off frequency (-3 dB) 1.4 GHz (±15%)

Time Response

Rise/fall time (10 % - 90 %)

250 ps (±15%)

Input

Noise equivalent power (NEP)

32 pW/\sqrt{Hz} (@ 760 nm, 100 MHz)

370 µW AC (for linear amplification, @ 760 nm) Optical saturation power

Detector

Detector Active area (FS/FST version)

10 mW CW (to prevent saturation, @ 760 nm)

Active area (FC version)

Si-PIN photodiode Ø 400 µm integrated ball lens

suitable for fibers up to 400 µm core diameter

Spectral range Max. sensitivity 320 - 1000 nm 0.51 A/W typ. (@ 760 nm)

Output

Output voltage range 1.9 V_{PP} (@ 50Ω output load)

for linear operation and low harmonic distortion

2.5:1 (@ f < 2.5 GHz) Output VSWR Output return loss 7.3 dB (@ f < 2.5 GHz) 50 Ω (terminate with 50 Ω load) Output impedance

Output noise 3.6 mV_{RMS} (24 mV_{PP}) typ. (@ 50 Ω load, no signal on

detector, measurement bandwidth 4 GHz)

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

HSA-X-S-1G4-SI_R11/TH,JMa/31MAR2022 Page 3 of 7


Ultra High Speed Photoreceiver with Si-PIN Photodiode

Specifications (continued)		
Input Flange	Material	1.4305 stainless steel, nickel-plated (FST flange) AIMg4.5Mn, nickel-plated (FS flange)
Coupler Ring (FST version only)	Material	1.4305 stainless steel, glass bead blasted
Power Supply	Supply voltage Supply current	+15 V 130 mA (depends on operating conditions, recommended power supply capability min. 200 mA)
Case	Weight	133 g (0.29 lbs) HSA-X-S-1G4-SI-FST incl. coupler ring 120 g (0.26 lbs) HSA-X-S-1G4-SI-FS 110 g (0.24 lbs) HSA-X-S-1G4-SI-FC
	Material	AlMg4.5Mn, nickel-plated
Temperature Range	Storage temperature Operating temperature	−30 °C +85 °C 0 °C +60 °C
Absolute Maximum Ratings	Optical input power (CW) Power supply voltage	12 mW (averaged) 20 V
Connectors Input Output Power supply	HSA-X-S-1G4-SI-FST 1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories	
		HSA-X-S-1G4-SI-FS 25 mm dia. unthreaded flange for free space applications
		HSA-X-S-1G4-SI-FC FC fiber optic connector (fix/permanent, FC/PC and FC/APC compatible)
	Output	SMA jack (female)
	Power supply	LEMO® series 1S, 3-pin fixed socket (mating plug type: FFA.1S.303.CLAC52)
		PIN 2 PIN 1 PIN 1 PIN 1: +15 V PIN 2: NC PIN 3 GND PIN 3: GND
Scope of Delivery	HSA-X-S-1G4-SI, internally threaded coupler ring (FST version only), LEMO® 3-pin connector, datasheet, transport package	
Ordering Information HSA-X-S-1G4-SI-FST HSA-X-S-1G4-SI-FS HSA-X-S-1G4-SI-FC	1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories.	
	HSA-X-S-1G4-SI-FS	25 mm dia. unthreaded flange for free space applications.
	HSA-X-S-1G4-SI-FC	FC fiber optic connector (fix/permanent, FC/PC and FC/APC compatible).

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E M T O

Datasheet HSA-X-S-1G4-SI **Ultra High Speed Photoreceiver** with Si-PIN Photodiode Spectral Responsivity 0.6 0.5 Photo Sensitivity (A/W) 0.4 0.3 0.2 0.1 0 200 300 400 500 600 700 800 900 1000 1100 Wavelength (nm) DB-Sens-HSA-X-S-1G4-SI_R01 HSA-X-S-1G4-SI-FST (1.035"-40 threaded free space input) **Dimensions** 43



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY


F E T O

DZ-HSA-X-S-1G4-SI_FST_R1

Ultra High Speed Photoreceiver with Si-PIN Photodiode

Dimensions (continued)

HSA-X-S-1G4-SI-FS (25 mm dia. unthreaded free space input)



DZ-HSA-X-S-1G4-SI_FS_R1

all dimensions in mm unless otherwise noted

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E T O




Datasheet HSA-X-S-1G4-SI

Ultra High Speed Photoreceiver with Si-PIN Photodiode

Dimensions (continued)

HSA-X-S-1G4-SI-FC (FC fiber optic connector)



DZ-HSA-X-S-1G4-SI FC R1

all dimensions in mm unless otherwise noted

FEMTO Messtechnik GmbH Klosterstr. 64 10179 Berlin · Germany Phone: +49 30 280 4711-0 Fax: +49 30 280 4711-11 Email: info@femto.de www.femto.de Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

© by FEMTO Messtechnik GmbH · Printed in Germany

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E M T O

HSA-X-S-1G4-SI_R11/TH,JMa/31MAR2022 Page 7 of 7